

JOBST® Confidence Measurement Form

Custom Made Lower Limb form



Patient: _____
 Surname First Name

Male Female

Date: _____

Compression Class	CL 1	CL 2	CL 3
Left			
Right			

Style


- AD B1T
 AG ET
 AG-HT AB1
 AT CT
 BT

Colour


- Beige
 Black
 Caramel
 Red Heather
 Anthracite Heather
 Jeans Heather

Additional Features

- Soft Fit (AD only)
 Silicone dots (AG only)
 Slip Form
 Comfort Zone Ankle
 Comfort Zone Knee
 Straight Foot
 Oblique Foot
 Closed Toe
 Open Toe
 Hallus Valgus Zone
 Left foot Right foot
 Deco Line
 Initials _____ (max 2 letters)
 Male Fly
 5cm Flowered Waist Band



Circumferences



Lengths

- K3 _____
 Measured from the waist (T), to the middle of the crotch without tightening.
 TT _____
 Measure from the waist (T), via the crotch to the back of where the garment will end. Do not apply tension.
 /T _____
 /H _____
 /K _____
 Approx. 2cm above the crotch/on the pubic bone.
 /G _____
 /F _____
 /E' _____
 /E _____
 /D _____
 /C _____
 /B' _____
 /B _____

Description of how hard to pull the tape measure

Tension	Description
0	Skin tension
+	Light tension
++	Firm tension
+++	Max tension

TIPS: cD and cG can be slightly firmer than Elvarex

Length measurements:

straight measurements on the inside of the leg.

Length /H:

straight measurements from the floor to the widest hip.

Length /T:

straight measurements from the floor to the end of the garment (straight measure to G and then contour G-H and H-T)

Foot lengths

- Length measurements in standing
- Circumferences in lying

Foot lengths

- _____ Total foot length
 _____ Medial length
 _____ Lateral length

